

## Solutions:

- 1,5% CTAB: Tris (75 mM), Cetyltrimethylammonium bromide (15 g/l) , EDTA (15 mM, pH8.0), NaCl (1.05 M), pH 8.0
- Chloroform:isoamylalcohol (24:1)
- **10xTE**  
100 mM Tris-Cl (desired pH)  
10 mM EDTA (pH 8.0)  
Sterilize solutions by autoclaving for 20 minutes at 15 psi (1.05 kg/cm<sup>2</sup>) on liquid cycle. Store the buffer at room temperature.
- **Tris-Cl**  
Dissolve 121.1 g of Tris base in 800 ml of H<sub>2</sub>O. Adjust the pH to the desired value by adding concentrated HCl.  
*pH HCl*  
7.4 70 ml  
7.6 60 ml  
**8.0 42 ml**  
(1 M) Allow the solution to cool to room temperature before making final adjustments to the pH. Adjust the volume of the solution to 1 liter with H<sub>2</sub>O. Dispense into aliquots and sterilize by autoclaving.  
If the 1 M solution has a yellow color, discard it and obtain Tris of better quality. The pH of Tris solutions is temperature-dependent and decreases approx. 0.03 pH units for each 1°C increase in temperature. For example, a 0.05 M solution has pH values of 9.5, 8.9, and 8.6 at 5°C, 25°C, and 37°C, respectively.
- **EDTA** (ethylenediaminetetraacetic acid)  
To prepare EDTA at 0.5 M (pH 8.0): Add 186.1 g of disodium EDTA•2H<sub>2</sub>O to 800 mL of H<sub>2</sub>O. Stir vigorously on a magnetic stirrer. Adjust the pH to 8.0 with NaOH (~20 g of NaOH pellets). Dispense into aliquots and sterilize by autoclaving. The disodium salt of EDTA will not go into solution until the pH of the solution is adjusted to ~8.0 by the addition of NaOH.

**Procedure:**

1. Cut part of a leaf (~2 cm) into small pieces, freeze in liquid nitrogen and grind leaf tissue to powder in a 2 ml safe-lock reaction tube containing a metal bead with a TissueLyser (Qiagen, Hilden, Germany)
2. Add 900  $\mu$ l boiled 1.5% CTAB buffer and incubate for 1h at 65°C.
3. Add 630  $\mu$ l of Chloroform:isoamylalcohol solution (24:1), and shake horizontally for 15 min at 75 rpm at room temperature.
4. Centrifuge at 15,000 rpm in a table centrifuge for 10 minutes.
5. Transfer the supernatant containing the nucleic acids into a new 2 ml reaction tube and precipitate DNA with 2/3 v/v ice-cold isopropanol (e.g. supernatant 700  $\mu$ l, isopropanol 450  $\mu$ l).
6. Shake gently and centrifuge at 15,000 for additional 10 min. Remove supernatant (be careful to keep the pellet).
7. Wash the precipitated nucleic acids with 1 ml 70% ethanol, centrifuge and remove supernatant (be careful to keep the pellet).
8. Dry in a vacuum centrifuge for 15 min.
9. Dissolve DNA in 50  $\mu$ l 1/10 TE buffer containing 5  $\mu$ g RNase A for RNA digestion. Keep the sample 2 hours at room temperature before freezing at -20°C.