

cDNA Synthesis Protocol NEB MuLV-Reverse transcriptase

1- Step 1:

1. Mix per reaction:
 - a. 1 μ l dNTP (10 mM)
 - b. 0.4 μ l Oligo (dT) 500 ng/ μ l (=100 μ M)
2. Transfer 1.4 μ l into each PCR tube.



3. Add 1 μ g RNA (volume to be calculated based on concentration of each sample) and adjust total volume to 16 μ l with RNase-free H₂O
4. 65 °C 5 min (temperature can be increased to max. 80 °C)
5. Put samples on ice, temperature of PCR block can be adjusted to 42 °C.



2- Step 2:

1. adding 4 μ l of the following master mix to samples kept on ice:

Reverse transcriptase buffer:	2.00 μ l
Reverse transcriptase (MULV):	0.25 μ l
RNase inhibitor:	0.50 μ l
RNA free water:	1.25 μ l

2. Continue PCR:
42 °C 60 min
90 °C 10 min
12 °C forever