

Chlorophyll and Carotenoid Determination (after Lichtenthaler 1987), a practical instruction.

Dear colleagues, the determination method for the photosynthetic leaf pigments, chlorophylls *a* and *b* and total carotenoids $x+c$, had been improved with freshly isolated pure chlorophylls and pure carotenoids. Their absorption spectra were measured in purest organic solvents using modern two-wavelength spectrophotometers which allowed to re-determine their specific absorption coefficients. The re-determined new equations for the quantitative determination of the photosynthetic pigments are found in the book series of *Methods in Enzymology*. Perhaps you have access to this book which is available in most international libraries. This is the method that most plant scientists apply for the quantitative determination of chlorophylls *a* and *b* in plant leaf extracts obtained with organic solvents such as acetone or methanol. This method also allows to simultaneously determine the level of chlorophylls *a* and *b* and that of total carotenoids in the same leaf extract solution.

LICHTENTHALER, H.K., Chlorophylls and carotenoids, the pigments of photosynthetic biomembranes. In: Douce, R. and Packer, L. (eds.), *Methods Enzymol.* 148, 350-382, Academic Press Inc., New York 1987.

You best extract your leaves or other green plant tissue with 100 % acetone. Centrifuge your extract at room temperature for 3 to 5 min at 300 to 500 x g so it will be **fully clear**. The latter is an absolute requirement, otherwise you will obtain wrong pigment values. Then you can determine the Chl *a+b* and total carotenoids ($x+c$) content in a two-wavelength spectrophotometer by using the equations given below. You measure the optical density, the absorption *A*, at 661.6 nm, at 644.8nm and at 470 nm, and then calculate with the equations the pigment amounts. These are indicated in μg per ml pigment extract solution. Then you need to consider the total extract solution you made. If you have 5 ml or 10 ml total extract solution, then you need to multiply by 5 or 10, respectively.

For 100 % pure acetone the new equations are given here. For other solvents check the original literature reference in *Methods in Enzymology* 1987 and cite this paper when you publish.

Chlorophyll *a*: $C_a = 12.25 A_{663.2} - 2.79 A_{646.8}$ (μg per ml solution)

Chlorophyll *b*: $C_b = 21.50 A_{646.8} - 5.10 A_{663.2}$ (μg per ml solution)

Total carotenoids: $C_{x+c} = (1000 A_{470} - 1.82 C_a - 85.02 C_b) / 198$ (μg per ml solution)

If your spectrophotometer does not allow to set exactly to 661.6 nm use 662 nm and for 644.8 nm use 645 nm. Best results are obtained when you measure at absorbance readings between 0.3 and 1.0 for the red chlorophyll peaks.

You will find the same equations with more detailed instructions and with examples for leaf chlorophyll and carotenoid levels in:

LICHTENTHALER, H.K. and C. BUSCHMANN: Chlorophylls and carotenoids – Measurement and characterisation by UV-VIS. *Current Protocols in Food Analytical Chemistry (CPFA)*, (Supplement 1), F4.3.1 - F 4.3.8 (2001) (John Wiley, New York).

This paper gives in Table F4.33 examples of the Chl and carotenoid content in leaves, as well as the weight ratios Chl *a/b* and Chl *a+b* /carotenoids, also expressed as $(a+b)/(x+c)$. Please

determine the pigment content on a leaf area unit (mg Chl a+b per m² or as µg cm⁻²) and also on a mg per 1g dry weight basis. Do avoid fresh weight as a reference system, because the water content of your leaf can change during the experiment. Therefore, fresh weight is not a valid reference system!

With the equations given in the 1987 *Methods in Enzymology* paper you can determine the levels of Chl a and b, as well those of total carotenoids (x+c) in the same pigment extract solution. 100% acetone is a very suitable solvent system for leaves. For isolated chloroplasts you better use 80 % acetone. Make sure that your extract solutions are fully clear, the best is to centrifuge the pigment extracts for 3 to 5 min at 300 x or 500 x g ! Pay attention: Any turbidity and light scattering in your extract solution yields wrong pigment values, in particular too low Chl a/b ratios are obtained, because the Chl b content in turbid solutions is estimated much too high. And the total carotenoid x+c levels are also fully wrong!

A further paper gives you information on the ways of extraction of Chls and carotenoids:

LICHTENTHALER, H.K. and C. BUSCHMANN: Chlorophylls and Carotenoids - Extraction, Isolation and Purification. *Current Protocols in Food Analytical Chemistry (CPFA)*, (Supplement 1), Unit F4.2.1-F4.2.6 (2001) (John Wiley, New York).

Examples for typical chlorophyll and carotenoid levels in leaves:

These are found in the following three papers:

Lichtenthaler H.K. Biosynthesis, accumulation and emission of carotenoids, α-tocopherol, plastoquinone and isoprene in leaves under high photosynthetic irradiance. *Photosynth. Research* 92: 163-179 (2007).

Sarijeva G., Knapp M., Lichtenthaler H.K. Differences in photosynthetic activity, chlorophyll and carotenoid levels, and in chlorophyll fluorescence parameters in green sun and shade leaves of *Ginkgo* and *Fagus*. *J. Plant Physiology* 164, 950 - 955 (2007).

Schindler, C., Reith, P. and Lichtenthaler, H.K.: Differential levels of carotenoids and decrease of zeaxanthin cycle performance during leaf development in a green and an aurea variety of tobacco. *J. Plant Physiol.*143, 500-507 (1994).

See and check chlorophyll and carotenoid levels also in these papers:

BABANI, F. and H. K. LICHTENTHALER: Light-induced and age-dependent development of chloroplasts in etiolated barley leaves as visualized by determination of photosynthetic pigments, CO₂ assimilation rates and different kinds of chlorophyll fluorescence ratios. *J. Plant Physiol.* 148, 555-566 (1996).

SCHINDLER, C. and H. K. LICHTENTHALER: Photosynthetic CO₂ assimilation, chlorophyll fluorescence and zeaxanthin accumulation in field-grown maple trees in the course of a sunny and a cloudy day. *J. Plant Physiol.* 148, 399-412 (1996).

LICHTENTHALER, H.K. and F. BABANI: Light adaptation and senescence of the photosynthetic apparatus: changes in pigment composition, chlorophyll fluorescence parameters and photosynthetic activity during light adaptation and senescence of leaves. In: Papageorgiou G. and Govindjee (eds.), *Chlorophyll Fluorescence: A Signature of Photosynthesis* (Chapter 30), pp.713- 736. Springer, Dordrecht, 2004.

Good luck with your measurements!

Karlsruhe, Germany, October 2010. Hartmut Lichtenthaler

Prof. Dr. Hartmut Lichtenthaler, Botany II, Karlsruhe Institute of Technology (KIT), University Division, Kaiserstr. 12, D-76133 Karlsruhe, Germany; Email: hartmut.lichtenthaler@kit.edu